

Crysta'Days | November 2023

Upscaling with the SCT-PLOT unit

**Secoya**  
FLUIDIFY PHARMA



A large, light teal shape with a rounded right side and a thin teal border occupies the upper and middle portions of the slide. In the bottom right corner, there are several overlapping teal shapes, including a triangle and a rounded rectangle, all in a darker teal color.

# 1

The new SCT-PLOT



# Secoya Crystallization Technology

SCT-pilot instrument

## SCT-PLOT uses:

- **Identical parameter settings** as the SCT-LAB environment
- **Single use fluidics** for easy and complete cleaning
- **All inserts and reactors are interchangeable between SCT-LAB and SCT-PLOT**
- **In-depth testing** continuous operation and **erroneous operation control**
- **Recipe determination**
- **GLP** version available
- Government funding for independent cGMP development

## Enhanced study delivers:

- **Deal** coupling with maturing and isolation strategy
- Full recipe for operations with industrial equipment

**WITH LIMITED USE of TIME and material**





# Secoya Crystallization Technology

SCT-pilot instrument

## SCT-PLOT changes in between versions:

- Selection of new pump type to increase cleaning capacity and reduce material losses: 150 mL dead volume reduced to 5 mL
- Complete emptying of the vessels
- Materials in contact with chemicals:
  - Glass and PTFE
  - Interior part of pumps in stainless steel metal, other materials on request
- Identical and validated flow rate stability






# Secoya Crystallization Technology

SCT-pilot instrument

## SCT-PLOT changes in bewteen versions:

- Selection of new pump type to increase cleaning capacity and reduce material losses: 150 mL dead volume reduced to 5 mL
- Complete emptying of the vessels
- Materials in contact with chemicals:
  - Glass and PTFE
  - Interior part of pumps in stainless steel metal, other materials on request
- Identical and validated flow rate stability






# Secoya Crystallization Technology

SCT-pilot instrument

## SCT-PLOT changes in between versions:

- Selection of new pump type to increase cleaning capacity and reduce material losses: 150 mL dead volume reduced to 5 mL
- Complete emptying of the vessels
- Materials in contact with chemicals:
  - Glass and PTFE
  - Interior part of pumps in stainless steel metal, other materials on request
- Identical and validated flow rate stability






# Secoya Crystallization Technology

SCT-pilot instrument

## SCT-PLOT technical specifications:

- Three different temperature zones from 5 to 85°C
- Flow rates between 1 and 100 mL/min, different pump heads may lead to better stability in different ranges
- Pressure drop possible up to 40 bar – very important for highly viscous materials
- Single use consumables
- Identical set of inserts and reactors
- Standard equipped with 2L glass lined vessels, agitated – may be exchanged with other vessels if required
  
- Pressure control for blockages
- Known OEM suppliers: Hubert thermostats and Knauer pumps.



# 2

## Case studies








# Secoya Crystallization Technology


SCT-pilot instrument long-term operation

## Adipic acid



Laboratory conditions: searching for different crystal attributes

Parameter	Adipic acid
Insert	Cooling - 0
Reactor	3 mL
Solution temperature (°C)	60
Solution concentration (mg/mL)	60
Flow rate (mL/min)	30 and 40
Nucleation temperature °C	various
Crystal growth conditions	Gathering at RT, 15 minutes maturing




remark: lower nucleation temperatures resulted in fast blockages during cleaning



# Secoya Crystallization Technology


SCT-pilot instrument long-term operation

## Adipic acid



Laboratory conditions: looking at crystal growth rate

Parameter	Adipic acid
Insert	Cooling - 0
Reactor	3 mL
Solution temperature (°C)	60
Solution concentration (mg/mL)	60
Flow rate (mL/min)	30
Nucleation temperature °C	20
Crystal growth conditions	Filtration at different collection times




remark: product gathered and filtered after one night at RT results in identical size



# Secoya Crystallization Technology

SCT-pilot instrument long-term operation

## Adipic acid



Laboratory conditions duplicated on pilot unit:  
Preparation of 8 L solution

Parameter	Adipic acid
Insert	Cooling – 0
Reactor	3 mL
Solution temperature (°C)	60
Solution concentration (mg/mL)	60
Flow rate (mL/min)	30
Nucleation temperature °C	17.5
Crystal growth conditions	Filtration 2 minutes after gathering

- Sampling was done every 30 minutes during 10 minutes to collect 300 mL of slurry,
- Slurry was gathered in bottle and filtered after 120 sec





# Secoya Crystallization Technology

## SCT-pilot instrument long-term operation

Laboratory conditions duplicated on pilot unit:  
Preparation of 8 L solution

Parameter	Adipic acid
Insert	Cooling - 0
Reactor	3 mL
Solution temperature (°C)	60
Solution concentration (mg/mL)	60
Flow rate (mL/min)	30
Nucleation temperature °C	17.5
Crystal growth conditions	Filtration 2 minutes after gathering




**4 hour run without disruptions**  
**Symetry in graph indicates stability of product over time**  
**Lab parameters fully extrapolated**



# Secoya Crystallization Technology


SCT-pilot instrument long-term operation

## Lactose



### Laboratory conditions set for **lactose**

Parameter	Lactose
Insert	Specific execution for highly viscous materials
Reactor	7 mL
Solution temperature (°C)	85
Solution concentration (mg/mL)	<b>850</b>
Flow rate (mL/min)	20
Nucleation temperature °C	5
Crystal growth conditions	Gathering in waiting tank for 18 hrs




remark: The insert and reactor are heavily modified to counter the viscosity effect



# Secoya Crystallization Technology


SCT-pilot instrument long-term operation

## Lactose



### Laboratory conditions set for **lactose**

Parameter	Lactose
Insert	Specific execution for highly viscous materials
Reactor	7 mL
Solution temperature (°C)	85
Solution concentration (mg/mL)	<b>850</b>
Flow rate (mL/min)	20
Nucleation temperature °C	5
Crystal growth conditions	Gathering in waiting tank for 18 hrs




remark: The insert and reactor are heavily modified to counter the viscosity effect



# Secoya Crystallization Technology

## SCT-pilot instrument long-term operation


### Lactose




Laboratory conditions set for **lactose**

Parameter	Lactose
Insert	Specific execution for highly viscous materials
Reactor	7 mL
Solution temperature (°C)	85
Solution concentration (mg/mL)	<b>850</b>
Flow rate (mL/min)	20
Nucleation temperature °C	5
Crystal growth conditions	Gathering in waiting tank for 18 hrs

remark: The insert and reactor are heavily modified to counter the viscosity effect






# Secoya Crystallization Technology

## SCT-pilot instrument long-term operation

Laboratory conditions set for **lactose**

Parameter	Lactose
Insert	Specific execution for highly viscous materials
Reactor	7 mL
Solution temperature (°C)	85
Solution concentration (mg/mL)	<b>850</b>
Flow rate (mL/min)	20
Nucleation temperature °C	5
Crystal growth conditions	Gathering in waiting tank for 18 hrs

remark: The insert and reactor are heavily modified to counter the viscosity effect









# Secoya Crystallization Technology

## SCT-pilot instrument long-term operation

### Laboratory conditions set for **lactose**


Parameter	Lactose
Insert	Specific execution for highly viscous materials
Reactor	7 mL
Solution temperature (°C)	85
Solution concentration (mg/mL)	850
Flow rate (mL/min)	20
Nucleation temperature °C	5
Crystal growth conditions	Gathering in waiting tank for 18 hrs

### Lactose



Axis = size in µm

### Lab repetition 1






# Secoya Crystallization Technology

## SCT-pilot instrument long-term operation


### Lactose



Laboratory conditions set for **lactose**


Parameter	Lactose
Insert	Specific execution for highly viscous materials
Reactor	7 mL
Solution temperature (°C)	85
Solution concentration (mg/mL)	850
Flow rate (mL/min)	20
Nucleation temperature °C	5
Crystal growth conditions	Gathering in waiting tank for 18 hrs

**2 hour run without disruptions**  
**Symmetry in graph indicates stability of product over time**  
**Lab parameters fully extrapolated**



Axis = size in µm

Pilot repetition 1





# Secoya Crystallization Technology



## SCT-pilot instrument long-term operation

### Laboratory conditions set for **lactose**


Parameter	Lactose
Insert	Specific execution for highly viscous materials
Reactor	7 mL
Solution temperature (°C)	85
Solution concentration (mg/mL)	850
Flow rate (mL/min)	20
Nucleation temperature °C	5
Crystal growth conditions	Gathering in waiting tank for 18 hrs

**Batch cross-test without passage through reactor results in 50% larger crystals with huge lumps of crystals**


### Lactose



Pilot repetition 1



Batch test





# Secoya Crystallization Technology



## SCT-pilot instrument long-term operation

Laboratory conditions set for **lactose**

Parameter	Lactose
Insert	Specific execution for highly viscous materials
Reactor	7 mL
Solution temperature (°C)	85
Solution concentration (mg/mL)	850
Flow rate (mL/min)	20
Nucleation temperature °C	5
Crystal growth conditions	Gathering in waiting tank for 18 hrs

**Batch cross-test without passage through reactor results in 50% larger crystals with huge lumps of crystals**

### Lactose



1 µm

Pil

100 µm

100 µm



# Secoya Crystallization Technology



## SCT-pilot instrument long-term operation

### Laboratory conditions set for **lactose**


Parameter	Lactose
Insert	Specific execution for highly viscous materials
Reactor	7 mL
Solution temperature (°C)	85
Solution concentration (mg/mL)	850
Flow rate (mL/min)	20
Nucleation temperature °C	5
Crystal growth conditions	Gathering in waiting tank for 18 hrs

**Batch cross-test without passage through reactor results in 50% larger crystals with huge lumps of crystals**


### Lactose



Pilot repetition 1



Batch test






# Secoya Crystallization Technology


## SCT-pilot instrument long-term operation

Paracetamol



Laboratory conditions set for **Paracetamol, directly applied to pilot**

Parameter	Paracetamol
Insert	T cross
Reactor	1m
Solution temperature (°C)	50
Solution concentration (mg/mL)	75
Flow rate (mL/min)	40/40
Nucleation temperature °C	RT
Crystal growth conditions	Collection in tank at 5°C for 30 min






# Secoya Crystallization Technology


SCT-pilot instrument long-term operation

## Paracetamol



Laboratory conditions set for **Paracetamol, directly applied to pilot**

Parameter	Paracetamol
Insert	T cross
Reactor	1m
Solution temperature (°C)	50
Solution concentration (mg/mL)	75
Flow rate (mL/min)	40/40
Nucleation temperature °C	RT
Crystal growth conditions	Collection in tank at 5°C for 30 min



A large teal number '3' is positioned on the left side of the slide. To its right, the text 'Production strategy' is written in a smaller, grey font. In the bottom right corner, there are several overlapping teal geometric shapes, including a triangle and a rounded rectangle, which serve as a decorative element.

# 3

**Production strategy**





# Secoya Crystallization Technology

SCT-ICE instrument

## INDUSTRIAL CRYSTALLIZATION UNIT – ICE is based on:

- **Identical Reactor setup** of laboratory equipment
- **cGMP design and execution**
- **Molecule specific or multiproduct approach**
- Reactors in parallel in standard equipment up to 10 modules
- **Cleaning in place**
- **Automated software, open for communication to its surroundings**
- Depending on solubilities and applied rates: 1.5 to 3 ton of solid product per module per year





# Secoya Crystallization Technology

SCT-ICE instrument

## INDUSTRIAL CRYSTALLIZATION UNIT – ICE is based on:

- **Identical Reactor setup** of laboratory equipment
- **cGMP design and execution**
- **Molecule specific or multiproduct approach**
- Reactors in parallel in standard equipment up to 10 modules
- **Cleaning in place**
- **Automated software, open for communication to its surroundings**
- Depending on solubilities and applied rates: 1.5 to 3 ton of solid product per module per year






# Secoya Crystallization Technology

Reliable Upscaling with identical set of parameters

Quality by Design approach (QbD):

Upscaling means placing **identical** reactors in **parallel**

Product overall quality is defined by temperature and flow stability






# Secoya Crystallization Technology

Reliable Upscaling with identical set of parameters


Quality by Design approach (QbD):

Upscaling means placing **identical** reactors in **parallel**

Product overall quality is defined by temperature and flow stability



Secondly, the reactors are opened







# Secoya Crystallization Technology

Reliable Upscaling with identical set of parameters

Quality by Design approach (QbD):

Upscaling means placing **identical** reactors in **parallel**

Product overall quality is defined by temperature and flow stability



Aim for the smallest possible operational zone

Real-time parameter control for fast adjustments and closing of reactors to maintain product quality





# Secoya Crystallization Technology

Reliable Upscaling with identical set of parameters

Quality by Design approach (QbD):

Upscaling means placing **identical** reactors in **parallel**

Product overall quality is defined by temperature and flow stability



Real-time parameter control for fast adjustments and closing of reactors to maintain product quality



# Secoya Crystallization Technology

Importance of operational zonewidth

Aim:


**Create identical residence time at identical temperatures**

Benefit:

- **More uniform conditions** of continuous operations means more **quality control by the design** of the instrument and operation
- Determination of **small operational zonewidths, extra space for adjustment of tests before running out of hand**

Test Flow Rate variation between  
**19.95 and 20.05 mL/min**

Boundary conditions:  
Quality of product is affected  
significantly outside of  
**19 and 22 mL/min**





# Secoya Crystallization Technology

Importance of operational zonewidth

Aim:


**Create identical residence time at identical temperatures**

Benefit:

- **More uniform conditions** of continuous operations means more **quality control by the design** of the instrument and operation
- Determination of **small operational zonewidths**, **extra space for adjustment of tests before running out of hand**

Test Flow Rate variation between  
**19.95 and 20.05 mL/min**

Boundary conditions:  
Quality of product is affected  
significantly outside of  
**19 and 22 mL/min**







# Secoya Crystallization Technology

Importance of operational zonewidth

Aim:


Create identical residence time at identical temperatures

Benefit:

- **More uniform conditions** of continuous operations means more **quality control by the design** of the instrument and operation
- Determination of **small operational zonewidths, extra space for adjustment of tests before running out of hand**

Test Flow Rate variation between **19.95 and 20.05 mL/min**

Boundary conditions:  
Quality of product is affected significantly outside of **19 and 22 mL/min**





# Secoya Crystallization Technology

Introduced into the parameter set on the CE equipment

Solution temperature

Peripheral system

Total pump flow rate

	Setpoint	Hystérésis	Stabilization time	Déviation time	Stabilization state	Remaining
<b>Thermostats</b>						
H1 :	65.0 °C	1.0 °C	600 sec	10 sec	<input type="radio"/> Hot tank temperature stabilized	0 sec
Tch2 :		1.0 °C	600 sec	10 sec	<input type="radio"/> Hot circuit temperature return stabilized	0 sec
H2 :	25.0 °C	10.0 °C	10 sec	10 sec	<input type="radio"/> Cold tank temperature stabilized	0 sec
Tf2 :		10.0 °C	10 sec	10 sec	<input type="radio"/> Cold circuit temperature return stabilized	0 sec
H3 :	10.0 °C	0.5 °C	600 sec	5 sec	<input type="radio"/> Cooling circuit of modules stabilized	0 sec
<b>Electrical heating</b>						
WH :	65.0 °C	1.0 °C	600 sec	10 sec	<input type="radio"/> Electrical heating temperature stabilized	0 sec
<b>Pressure valve</b>						
BPR CH :	9.0 bar	0.5 bar	300 sec	10 sec	<input type="radio"/> Hot circuit pressure stabilized	0 sec
BPR F :	12.0 bar	3.0 bar	10 sec	15 sec	<input type="radio"/> Cold circuit pressure stabilized	0 sec
<b>Pumps</b>						
P1 :	200.0 ml/min	1.0 ml/min	180 sec	5 sec		
P2 :	21000.0 ml/min	1.0 ml/min	15 sec	15 sec		

**Réservoir volume**    **Pipe volume**

V1 : 40.0 L    1.0 ml    **Reset Value**

V2 : 40.0 L    1.0 ml

**Setting to start sequence**



# Secoya Crystallization Technology

Introduced into the parameter set on the CCE equipment

Flow rate in reactor

Nucleation temperature

Inside modules

The screenshot displays the control interface for three modules (Module n°1, 2, and 3). Each module has a 'Setpoint flow' section and a 'Process value' section. The 'Setpoint flow hot circuit' for Module n°1 is set to 20.0 ml/min, with a hysteresis of 5.0 ml/min, a stabilization time of 17 sec, and a deviation time of 15 sec. The 'Process value' for Module n°1 shows a flow rate of 0.3 ml/min. The 'Setpoint flow cold circuit' for Module n°1 is also set to 20.0 ml/min, with a hysteresis of 5.0 ml/min, a stabilization time of 15 sec, and a deviation time of 15 sec. The 'Process value' for Module n°1 shows a flow rate of 0.2 ml/min. The 'Setpoint temperature module' for Module n°1 is set to 25.0 °C, with a hysteresis of 30.0 °C, a stabilization time of 15 sec, and a deviation time of 15 sec. The 'Process value' for Module n°1 shows a temperature of ##### °C. The interface includes a navigation bar at the bottom with buttons for 'Copy setpoint for all modules', 'Paramètres PID', 'Mod. 1', 'Mod. 2', 'Mod. 3', 'Mod. 4', 'Mod. 5', and 'Number of module : 5'. A status bar at the bottom shows 'Time' and 'Text' fields, along with icons for home, settings, stop, graph, and spray, and warning and power icons.




# Secoya Crystallization Technology

Reliable Upscaling with identical set of parameters

Upscaling means placing **identical** reactors in **parallel**

Product overall quality is defined by temperature and flow stability






# Secoya Crystallization Technology

Upscaling strategy – based on actual examples

- Production request: **25 kg/day** for a production of 100 kg/week
- Only working hours, operations possible during 12 hrs: 1 hr setup, 9 hrs to run, rest for cleaning (2 shift basis with overlap)
- Solubility in solvent: 250 mg/mL
- Solubility at equilibrium (5°C): 50 mg/mL
- Possible yield: 200 mg/mL solvent
- Obtained: 94 % yield = 196 mg per mL solution injected
- Flow rate: 30 mL solution/min/reactor
- Production rate per hour ~ **350 g/hr/reactor**
- 8 reactors produce per 9 hours 25.2 kg
- Use bonus of two modules to shorten production times




1 Growth tank to collect material at temperature, agitated  
Filtration with dimensionalized candle filter to shorten filtration time and to allow washing and during inside the tank with DrMueller



# Secoya Crystallization Technology

Upscaling strategy – based on actual examples

- Production request: **50 ton per year**, split in weekly productions
- 5 days' work week, cleaning procedure included, nucleation reactor may run overnight
- Solubility in solvent: 850 mg/mL
- Solubility at equilibrium (5°C): 150 mg/mL
- Possible yield: 700 mg/mL solvent
- Obtained: 95 % yield = 665 mg per mL solution injected
- Flow rate: 20 mL solution/min/reactor
- Production rate per hour ~ **798 g/hr/reactor**
- 50 reactors produce per 3 days 2.87 ton
- Ideal scenario: 17 batches per year which leaves space for events + increased production



3 growth tanks are filled consecutively to cope with 18 hr growth time  
Filtration and CIP + module cleaning on Thursday and Friday




# Secoya Crystallization Technology

Upscaling strategy – based on actual examples

- Production request: reactive crystallization of inorganic salt
- 5 days' work week, 2 shift basis for a 12 hr production + cleaning possibility
- **5 kg solids/hour** = 200 kg/week = 9 ton/45 weeks operation
- Possible yield: 100 mg/mL solution with 1/5 mixing ratio with antisolvent

!!!! This reactive crystallization does not require a reactor, only turbulent mixing and a larger execution of the insert is applied, running at 200 mL/min for the solution side and 1000 mL/min for the antisolvent side

- 4 nozzles applied, directly implemented onto a Nutsche filtration unit. Production rate attained **1.2 kg/hr/reactor**



# 4

**Other production strategies**  
**PIPELINE – client input**





## How to implement Secoya at industrial level - Discussion

We see levels of entry that are more feasible to introduce the technology into industry:

- HPAP as we control the size of product and manipulations are over
- Specific formulations: Drug product level, skipping one or more steps in DP development
- Benefit to couple to continuous tableting machines - Fette – at DP level



More classic AP substance production, splitted from Drug product manufacturing, we see interest, but difficulties in entering into these sites as the **change** necessary to implement our units is drastic/ 'adventorous'



# Preparation of high quality seeds – both dried or in suspension

## Seeding Strategy

- Generate seeds of desired size and PSD, at equilibrium condition
- ▣ Add seed suspension to solution – might be dried (storage) or in suspension
- ▣▣ Perform a 'classic' batch crystallization





# Produced slurry used as starting point for MSMPR setup

Use as continuous feeder of nuclei into MSMPR setup

- Generate seeds
- ▣ Continuously added to first CSTR tank
- ▣▣ Product is removed from first and sent to a second at lower temperature

